

We Serve, You Cure



AS Solutions Inc.



Plot No .1388, LIG I, 10th Street, TNHB Colony, Velachery, Chennai - 600042

Tel.fax: 044 - 22450321 Ph. No: +91 - 9566145898





# Solutions Inc.

# Your Respiratory Partner

# **NELLCOR OXIMAX N-560**

## Monitor with confidence.

The Nellcor™ OxiMax N-560 pulse oximeter from Covidien is a compact, affordable bedside monitor that delivers full-size performance. It incorporates the latest Nellcor™ signal processing technology for accurate, reliable SpO2 and pulse rate measurements even in conditions of low perfusion and signal interference, including patient motion. This technology, along with the full line of versatile OxiMax sensors, allows clinicians to confidently and effectively monitor a wide range of patients.



#### FEATURIN G LOSAT EXPANDED ACCURACY

The LoSat expanded accuracy feature of our OxiMax adhesive sensors assures you of a wide accuracy range (60% to 100% SpO2) when used with the Nellcor OxiMax N-560 pulse oximeter. This allows improved patient assessment at the more challenging lower saturation levels common in infants with congenital heart disease.

## NELLCOR™ OX IMAX N-560 PULSE OX IMETER FEATURES

- Monitors neonates through adults
- > Simple, intuitive operation and space-saving design.
- The Covidien distinctive digital signal processing technology enhances performance during low perfusion, patient motion and other forms of signal interference.
- Variable pitch beep tone enables clinicians to hear pointby-point changes in SpO2
- Display changes from green to red during alarms.
- Eight-hour battery life.
- 24-hour trend memory

hospital sleep screening.

- Optional paging system supplements bedside and central
- station alarms to reduce response time.
  Integrated Score analysis functionality stream lines in-

## NELLCOR™ OX IMAX SENSOR ADVANTAG ES

As part of the Covidien Nellcor<sup>TM</sup> OxiMax pulse oximetry system, the OxiMax N-560 pulse oximeter works with the full line of OxiMax sensors, which help expand your patient monitoring options.

- OxiMax SpO2 forehead sensor offers improved performance when patients have low perfusion.
- OxiMax nonadhesive sensors provide a gentle, singlepatient-use alternative to adhesive sensors for patients with fragile skin.

Tel.No: +91-44-22450321 Ph. No: +91 - 9566145898

Email: sales@amssolutionsinc.com

Plot No .1388, LIG I, 10th Street,

TNHB Colony, Velachery,





# Solutions Inc.

# Your Respiratory Partner

## **NELLCOR OXIMAX N-560**

### **PERFORMANCE**

### Measurement range

SpO2: 1% to 100%

Pulse rate: 20 to 250 beats per minute (bpm)

Pulse amplitude: 0.03% to 20%

#### Accuracy

Saturation (% SpO2 ± 1 SD)

Adult: 70% to 100% ± 2 digits

60% to 80% ± 3 digits

Neonate: 70% to 100% ± 2 digits

60% to 80% ± 4 digits

Low perfusion: 70% to 100% ± 2 digits

Pulse rate 20 to 250 bpm ± 3 digits

Low perfusion: 20 to 250 bpm ± 3 digits

ELECTRICAL

#### Instrument

Power requirements: 100 to 240 VAC, 25 to 39 VA
Fuse rating: Qty 2, 2A, 250 volts, slow-blow

IEC (5x20 mm)

Battery

Type: NiMH

Battery capacity: Minimum of 8 hours using new, fully

charged battery with no alarms, no

serial data

#### PHYS ICAL CHARACTERISTICS

Weight 3.07 lb (1.39 kg)

Size 2.95 H x 9.06 W x 5.04 D (inches)

(7.5 cm x 23.0 cm x 12.8 cm)

Equipment IEC 60601-1 IPX1

CSA 601.1 ISO 9919:2005

UL 2601-1 EN/IEC 60601-1-2 (2nd ed.)

## FEATURES & SPECIFICATIONS

### Emissions compliance

EN55011, Group 1, Class B

#### Output

- > EIA-232 digital output from serial port
- > OxiNet Communications Output Protocol
- > Nurse call capability

#### Display/indicators

- > High quality LED display changes from green to
- > red during alarm conditions.
- > Pulse amplitude indicator (10-segment LED).
- ➤ Visual indicators: Pulse search, audible alarms silenced or off, interference indicator, battery charging, data in sensor and Nellcor™ SatSec onds alarm management clock.

#### Alarms

- Audible and visual alarms for high/low satura tion and pulse rate.
- > SatSeconds Alarm Management settings: 10. 25. 50 and 100. or OFF

10, 25, 50 and 100, 61 OFF

Audible and visual warning indicators for low battery and sensor off.

Audible sensor disconnect alarm.

## Optional accessories

- ▶GCX mounting adapter
- > Pole mount clamp
- > Interface cables
- > Visual alarm indicator

Tel.No: +91-44-22450321 Ph. No: +91 - 9566145898

Email: sales@amssolutionsinc.com Plot No .1388, LIG I, 10th Street,

TNHB Colony, Velachery,

Chennai - 600 042.