

We Serve, You Cure



AS Solutions Inc.



Plot No .1388, LIG I, 10th Street, TNHB Colony, Velachery, Chennai - 600042

Tel.fax: 044 - 22450321 Ph. No: +91 - 9566145898





Your Respiratory Partner

NEW PORT HT50

MONITORED PARAMETERS

Paw (Peak, Mean, Base): VT (insp. tidal volume)
Internal Battery: f total (total frequency)

ALARMS

Paw(High Airway Pressure):4 to 99 cmH2O/4 to 99 mbar Paw(Low Airway Pressure):3 to 98 cmH2O/3 to 98 mbar

High Baseline Pressure: Automatically set
Low Baseline Pressure: Automatically set
V_I(High Insp. Minute Volume): 1.1 to 50 L
V_I(Low Insp. Minute Volume): 0.1 to 49.0 L

Apnea (Paw/Apnea): 30 sec
Device Alert: malfunction alert
Power Switchover: AC to DC

Low Battery:minimum 30 minutes of back up power

remains

Silence/Reset: 60 sec.

*PEEP + support cannot exceed 60 cmH2O/mbar

SPECIFICATIONS

MISCELLANEOUS

Operating Temperature: -18°C to 50°C
Operating Altitude:0 to 15,000 ft (0 to 4,572m)

Storage Temperature: -40°C to 65°C RS-232: 8 pin Semconn connector

19,200 Baud Remote Alarm/Nurse Call can be set up normally open or closed .

ORDERING INFORMATION

HT50-H1-B: Ventilator

HT460231: Compact stand with support arm
HT460700: Air/Oxygen Entrainment Mixer

HT600044: Oxygen Blending Bag Kit

HT460085: Automotive lighter 12V Power cord

HTECH01: E-cylinder holder for stand

Tel.No: +91-44-22450321 Ph. No: +91 - 9566145898

Email: sales@amssolutionsinc.com

Plot No .1388, LIG I, 10th Street,

TNHB Colony, Velachery,

Chennai - 600 042.





Solutions Inc.

Your Respiratory Partner

NEW PORT HT50

SPECIAL FEATURES



- > Volume or Pressure Control Ventilation
- Invasive or Non-Invasive Ventilation
- Pressure Support
- Back-up Ventilation
- > PEEP/CPAP, Built-in
- Comprehensive Alarms/Monitoring
- Dual Pac Internal Battery System, up to 10 hours of operation
- > Unique Gas Delivery System Requires No External Air Compressor
- > Optional Air/Oxygen Mixer for High Pressure 0.21-1.00
- Optional Blending Bag for Low Flow Systems
- Secondary back-up battery built-in for added safety

PATIENT RANGE

10 kg Infant, Pediatric and Adult

CONTROL SETTINGS

Modes (Pressure or Volume):A/CMV SIMV SPONT

Volume Control (Tidal volume): 100 to 2,200 mL

f(Frequency):1 to 99 b/min

tl(Inspiratory Time):1.1 to 3.0 sec

PEEP/CPAP: 0 to 30 cmH2O/0 to 30 mbar

Pressure Support: 0 to 60 cmH2O/0 to 60 mbar*

Flow: 6 to 100 L/min

I:E Ratio: 1:99 to 3:1

Pressure Control

(Target Pressure): 5 to 60 cmH2O/5 to 60 mbar

P trig (Sensitivity): -9.9 to 0 cmH2O/-9.9 to 0 mbar

Panel Lock: ON/OFF

Manual Inflation: 3 sec maximum

(single operator

controlled breath)

FIO2: .21 to 1.00

SPECIFICATIONS

PNEUMATICS

Gas Delivery System: internal dual micro-piston

Air/Oxygen Mixer (Optional):35 to 90 psig oxygen

POWER REQUIREMENTS :100-240 VAC, MAX. 2A

PNEUMATICS

:Primary lead acid, 12V 5AH when new and fully

charged, provides up to 10 hrs of operation

:Secondary NiMH, 12V 2AH provides emergency

back up power for minimum 30 mins.

DIMENSIONS

Width: 10.63 in (27 cm)

Depth: 7.87 in (20 cm)

Height: 10.24 in (26 cm)

Tel.No: +91-44-22450321 Ph. No: +91 - 9566145898

Email: sales@amssolutionsinc.com

Plot No .1388, LIG I, 10th Street,

TNHB Colony, Velachery,

Chennai - 600 042.